

Investigating the Oxidation of a Cobalt-based Catalyst Using X-ray Photoelectron Spectroscopy

Chris Baily, Paul Mack, Thermo Fisher Scientific, East Grinstead, UK

Key Words

K-Alpha, Cobalt Oxidation, Hydrodesulfurization Catalysis, Metal Oxides, Surface Analysis, XPS

The Thermo Scientific K-Alpha system was used to investigate the surface chemical change between a fresh γ -alumina supported Co(Ni)MoS hydrodesulfurization catalyst and the same catalyst after exposure to air. By using X-ray photoelectron spectroscopy (XPS) identification and quantification of chemical differences between the two samples was possible.

Introduction

Chemical state determination and quantification are of critical importance when understanding the performance of catalysts. Metal-based catalysts are typically used for redox reactions, and so careful control of the metal oxidation state is required to ensure that the catalyst performs successfully. The surface of the material is where the reaction occurs; thus a method of measuring the surface chemistry of the catalyst is required. X-ray photoelectron spectroscopy is an excellent technique for this type of characterization. Samples are analyzed under vacuum, and so surface changes during the analysis are kept to a minimum. XPS is an extremely surface sensitive technique too. Spectra result from within the top 10 nm of the surface, which allows XPS to focus on the crucial surface area. Finally, XPS can both identify and quantify both elemental compositions and the chemical states of the elements detected.

The high performance of the Thermo Scientific K-Alpha surface analysis system enables researchers to use XPS to quickly and clearly determine the composition and chemical states of the surfaces of catalyst samples. In this application note we describe how the K-Alpha was used to analyze a γ -alumina supported Co(Ni)MoS catalyst to obtain chemical state information from both a fresh and air exposed sample and determine the type and level of chemical change the catalyst underwent when exposed to air.



Method

The catalyst pellet was removed from its storage under hexane and mounted directly onto double-sided copper tape, which had been adhered to the sample platter. The sample block was then immediately loaded into the fast entry load-lock, and was automatically pumped down and transferred into the analysis chamber. The pellet was analyzed for elemental and chemical state information. Once this analysis was complete, the sample was transferred out of the system and exposed to air at atmospheric pressure for 3 minutes. The sample was then pumped down and transferred back into the analysis chamber. The subsequent XPS analysis was compared with the original results to determine the level of oxidation of the sample surface following exposure to air.

Results

The cobalt $2p_{3/2}$ spectra recorded on fresh and aged catalysts are shown in Figures 1 and 2 respectively. The different phases of Co were quantified using the integrated peak fitting software in the Avantage data system according to reference 1.

The spectrum was fitted using a multicomponent fit assuming Co_9S_8 , CoMoS and Co(II) oxide phases. No significant oxidation was observed when transferring the sample from container to K-Alpha, confirming that the sample had the expected composition.

More oxidation was observed after exposure to air, leading to an increase in the relative amounts of Co(II) oxide and Co_9S_8 at the surface and a decrease in the relative amount of CoMoS .




Figure 1: $\text{Co}2p_{3/2}$ spectrum of fresh catalyst




Figure 2: $\text{Co}2p_{3/2}$ spectrum of aged catalyst




Figure 3: Mo3d spectrum of fresh catalyst

A similar increase in oxidation is detected in the molybdenum spectra, shown in Figures 3 and 4.

As with the cobalt spectra, more oxidation was observed after exposure to air, leading to an increase in the relative amounts of Mo(IV) oxide and of Mo(VI) oxide at the surface and a decrease in the relative amount of Mo(II) sulfide.

Table 1 shows the comparison table of the full chemical state quantification for the fresh and aged samples.

The quantification shows that the amount of alumina support has stayed the same and that the sulphur, molybdenum and cobalt have all become more oxidized. The identification of increased levels of oxidation enables the analyst to determine the quality of the surface of the catalyst.

Summary

The K-Alpha XPS system can be used to quickly and accurately identify and quantify the chemical states of the elements present in catalyst samples.




Figure 4: Mo3d spectrum of aged catalyst

XPS Peak	Chemical State	Fresh Sample (At. %)	Aged Sample (At. %)
Al2p	Al ₂ O ₃	23.97	24.65
S2p _{3/2}	Sulfide	7.01	5.34
S2p _{3/2}	Oxysulfide	0.30	0.90
Mo3d _{5/2}	Mo(II) Sulfide	2.62	1.43
Mo3d _{5/2}	Mo(IV) oxide	0.14	0.86
Mo3d _{5/2}	Mo(VI) oxide	0.24	0.74
C1s	C-C/C-O/C=O	22.79	17.81
O1s	O-C/O-metal	40.77	46.14
Ni2p	Ni(II) sulfide	1.64	2.13
Co2p _{3/2}	Co ₉ S ₈	0.07	0.10
Co2p _{3/2}	CoMoS	0.33	0.22
Co2p _{3/2}	Co(II) oxide	0.13	0.20

Table 1: Full chemical state quantification for both the fresh and aged sample

Acknowledgements

Thermo Fisher Scientific would like to thank Dr. Pavel Nikulshin, Samara State Technical University, Russia, for supplying the catalyst sample used in this work.

References

- Gandubert et al., *Oil & Gas Science and Technology* – Rev. IFP, Vol. 62 2007, No. 1, pp. 79-89

www.thermoscientific.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +27 11 822 4120

Australia +61 3 9757 4300

Austria +43 1 333 50 34 0

Belgium +32 53 73 42 41

Canada +1 800 530 8447

China +86 10 8419 3588

Denmark +45 70 23 62 60

Europe-Other +43 1 333 50 34 0

Finland/Norway/Sweden

+46 8 556 468 00

France +33 1 60 92 48 00

Germany +49 6103 408 1014

India +91 22 6742 9434

Italy +39 02 950 591

Japan +81 45 453 9100

Latin America +1 561 688 8700

Middle East +43 1 333 50 34 0

Netherlands +31 76 579 55 55

New Zealand +64 9 980 6700

Russia/CIS +43 1 333 50 34 0

Spain +34 914 845 965

Switzerland +41 61 716 77 00

UK +44 1442 233555

USA +1 800 532 4752

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific